

Opti-Lube®

Lubricating Instrument Rinse Solution



- Ready To Use
- Anticorrosive
- Non-Toxic
- Non-Flammable

Opti•Lube®

Lubricating Instrument Rinse Solution

-Ready to Use-

Anticorrosive Non-flammable Non-Toxic

Keeps surgical instruments new!



Description: Use on stainless steel surgical and dental instruments after cleaning and rinsing as a final rinse prior to steam sterilization. Effectively rinses instrument surfaces by removing mineral deposits left by tap water, detergent residues and residual soil left after cleaning. Applies a very light surgical grade lubrication to instrument surfaces that is absorbed into the instrument surface during the steam sterilization cycle. Prevents staining, corrosion and discoloration of instruments. Prolongs instrument life and helps moving parts to glide smoothly.



Applications: Simple to use. Spray instruments using *Opti•Lube* Lubricating Rinse Solution. Leave surfaces wet. Package and autoclave.



Safety: *Opti•Lube* is non-toxic, non-flammable and non-caustic. Independent laboratory studies show *Opti•Lube* Lubricating Instrument Rinse Solution to be non-irritating. No special handling is required. *Opti•Lube* Lubricating Instrument Rinse Solution is packaged following strict FDA Good Manufacturing Practices.



PRODUCT

16 oz. Trigger Spray Bottle
One Gallon Plastic Pour Bottle

REORDER

OLR12-016
OLR04-128

PACKAGING

12 / Case
4 / case

Distributed by:

For technical data and product information
contact:



Rolling Meadows, IL 60008 USA
Telephone: 1-888-253-2536
(Outside USA: 1-847-454-0835)